

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Vibration Measurement Trainer

Product Code :
SCHOOL-MES880016



Description :

Vibration Measurement Trainer

Technical Specification :

FEATURES:

The whole system is self-contained with all accessories included.

The instrument is enclosed in a tapered fiberglass box and all necessary test points are brought out with schematic circuitry for understanding the working of the instrument.

SPECIFICATION:

Vibration Source: Vibration generator capacity 5 new tons.

Power Amplifier: 0 to 10V P-to-P variable by potentiometer.

50 to 1000Hz Variable by potentiometer.

Indicator 4½ digit 16 X 2 line LCD display.

Parameter Displayed Acceleration, Velocity, Displacement and frequency.

Vibration Sensor: Accelerometer (piezo electric crystal).

Power required for operation

230V, 50Hz, 1 phase

230V, 60Hz, 1 phase

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080