

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Gasoline Direct Injection Engine Test Bench

Product Code :
SCHOOL-THM830011



Description :

Gasoline Direct Injection Engine Test Bench

Technical Specification :

FEATURES:

Trainer can be equipped with gasoline direct fuel injection, electronic ignition, system, and wideband oxygen sensor.

Depending on specific application, trainer comes with specific and customized options.

SPECIFICATIONS:

Demonstration of systems to groups of students without the access limitations of a complete vehicle.

Provide your students with the latest in vehicle systems and diagnostic technologies.

Engine systems respond to inserted faults with real world symptoms, OEM DTC's, scan tool Data Stream changes and malfunction Indicator light operation and oscilloscope pattern changes.

Depending on specific application trainer can be equipped with the latest in engine performance management systems and may include :

Gasoline Direct Injection

Variable Valve Timing

CAN Bus

EVAP

Electronic Throttle Control

Wide-Band Oxygen sensor

Catalytic Converter

Electronic ignition

OEM A/C system

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080